

# BMC<sup>®</sup> RESmart™ APAP Machine with Humidifier

The RESmart™ Auto (APAP) Machine includes advanced sensor technology to automatically track your respiratory condition.

## Features include...

- Flexible tri-mode operation (CPAP, APAP and Titration)
- Powerful onboard memory with i-Code (No SD Card required)
- 0 ~ 60 minute ramp time
- Automatic airflow on / off
- Automatic mask leak detection and compensation
- Mask Off alarm
- Display with backlight and illuminated user buttons
- Convenient carry handle
- Easy review of the last night's sleep data
- Removable heated humidifier



## Standard Kit



Technical Specifications Summary	
Size	220 × 194 × 112 mm   313 × 194 × 112 mm with humidifier
Weight	<1.6Kg   <2.4Kg with humidifier
Power	100~240VAC, 50/60Hz
Operating Temperature   Humidity	50°C ~ 40°C   <80% Non-Condensing
Storage Temperature   Humidity	-20°C ~ 55°C   <93% Non-Condensing
Pressure Range	4~20 hPa (4~20cmH20)
Sound Pressure Level	<30dB at a working pressure of 10hPa
Maximum Flow Rate	>35 litres per minute
Pressure Display accuracy	+/- 0.5hPa
Approvals	IEC 60601-1 General Requirements for Safety of Medical Electrical Equipment
	ISO17510 Sleep Aponea Breathing Therapy
	Australian Therapeutic Goods Administration (TGA)

*This material is provided for informational and promotional purposes only and must not be considered as a substitute for professional medical advice, diagnosis or treatment. Any decisions made as a result of viewing this document, or visiting our website should be acted or relied upon only after consultation with your doctor and/or sleep clinic.*