

Small, light and easy to wear



N5B Micro Nasal CPAP Mask

The ultra-lightweight and minimal design means the mask is easy to assemble and clean.

The innovative design means the BMC N5B Nasal Mask has more movement and flexibility, reducing mask leaks to a minimum

Most importantly it is light and comfortable to wear. The BMC N5B Fir Pack includes all 3 sizes to help ensure a perfect fit.



Widened vision, improved seal

The design without forehead frame is suitable for patients suffering from nasal septal deviation.

The dynamic adaptive cushion design maximises patient comfort whilst minimising patient disturbance caused by mask leaks during sleep.




N5B Micro Nasal Mask



N5B Micro Nasal Mask

Minimal contact

Thanks to the minimal contact design, less pressure is applied to the face, while maintaining a stable seal on the face.



Platform-based frame

The mask frame is compatible with all three cushion sizes with a simple one-click lock.

Real comfort


The groove design on the cushion plus the ball and socket elbow provides a stable mask experience, even for the most restless of sleepers.



Crown-style headgear

The crown-style headgear provides the user with a more stable mask wearing experience.

Technical Specifications



Product Codes

Description	Size	Code
N5B Nasal Mask	Large	10020319
N5B Nasal Mask	Medium	10020320
N5B Nasal Mask	Small	10020321
N5B Nasal Mask	All Size Kit (L, M & S)	10020322

Description	Size	Code
Cushion Assembly	Large	60380242
Cushion Assembly	Medium	60380243
Cushion Assembly	Small	60380244

Description	Size	Code
Headgear Assembly	/	61030060